

Frequently Asked Questions & Answers

Q) Can the BiteStrip™ be reused?

A) No, the BiteStrip is a single-use, disposable device. After the study ends, the result is permanently etched on the display and cannot be erased. The BiteStrip can be a part of the patient's permanent record.

Q) What is the shelf-life of the BiteStrip™?

A) The shelf-life of the BiteStrip is two years from the date it was manufactured. The expiration date is clearly marked on the package.

Q) Why is a doctor's prescription required for the BiteStrip™?

A) Bruxism is a medical disorder and should be diagnosed and treated by a dentist. There are many variables in determining the existence and frequency of bruxism. A final diagnosis and treatment recommendations can only be made by a qualified professional who will use both the clinical evaluation and the BiteStrip study results.

Q) How long does the result image last on the display?

A) The result image is permanent. After the green sticker is removed, the display is readable for at least five years. If the green sticker remains in place for more than three months, corrosion of the display may occur.